

P50

Accessories

- 1 injection arm including injection tubing system
- 1 stand (with 2 height-adjustable hanging devices)
- 2 bags of powder to mix artificial blood
- 1 infusion bottle with screw cap, 250ml
- 1 rubber cork with hole
- 1 suspension net for infusion bottle
- 1 infusion system with air filter
- 1 plastic cup
- 1 disposable syringe 5ml
- 2 injection cannulas for peripheral blood sampling, 20 and 21 gauge (recommended cannula size)
- 2 replacement tubing systems with white sealing caps
- 4 white replacement sealing caps
- 1 can of talcum powder

A few advance remarks

- The skin of the 3B Scientific® injection arm is made of an elastic material which retracts after small punctures so that the used injection site cannot be seen. After prolonged usage and many punctures in the same area it may be necessary to replace the skin. Before this is done, the tubing system should be emptied and cleaned as described below.
- Besides the puncturable veins (filled with fluid), additional vein structures are shown in the skin. These veins cannot be punctured, but they are represented because the arm is an anatomical replica. The following veins can be punctured: basilic vein, cephalic vein, median cubital vein, dorsal venous rete of hand.
- The tubing system is subject to wear and tear and is easy to replace. 2 spare tubing systems are supplied with the injection arm.
- Please use only injection cannulas sized 20 or 21 gauge.
- From time to time or when you are planning not to use the injection arm for several days it is recommended to rinse the tubing system and the infusion system with water in order to remove any artificial blood residues. To do this, fill the artificial blood from the infusion bottle into a jar and clean the system following the procedure described in the section "Mounting the injection arm", using water instead of artificial blood.
- An alternative to artificial blood is water or standard infusion bottles with salt solution. In these cases cleaning is required less frequently.

Mounting the injection arm

- Attach the upper hanging device at the red mark on the stand and the lower one at a height of approx. 10 cm. (The infusion bottle will later be hung into the upper hanging device. You can also hang it above the red mark to increase the pressure of the artificial blood. However, with the bottle hanging higher, fluid might escape from the puncture holes).
- Hang the injection arm onto the lower hanging device using the ring mounted in the end of the arm.

Injection Arm

English

- Fill the contents of one bag of artificial blood powder into the infusion bottle and fill up with 250 ml of water. Close the bottle and shake well so that the powder dissolves completely.
- Now replace the screw cap of the bottle with the rubber cork.
- Pull the suspension net over the infusion bottle (beginning at the bottle neck).
- Close the flow control of the infusion system on the infusion tube (fig. 1)
- Now push the tip of the infusion system as far as possible through the hole in the rubber cork (fig. 2)
- Hang the infusion bottle into the upper hanging device of the stand (marked red).
- Fill the drip chamber of the infusion system by pressing together the chamber walls several times and open the ventilation flap (fig. 3)
- Remove the seal at the end of the infusion tube.
- Both ends of the tubing system are visible at the end of the injection arm. This tubing system simulates the veins through which the artificial blood will later be conducted. First remove the white sealing cap of the short tube and slip the end of the tube over the end piece of the infusion tube (fig. 4)
- Then remove the white sealing cap of the long tube and hang the end of the tube into the plastic cup.
- Now open the flow control on the infusion tube (fig. 5). As soon as the artificial blood reaches the cup close the flow control again and replace the white sealing cap onto the end of the tube.
- Now open the flow control of the injection tube again, so that withdrawn artificial blood can flow back in.
- You can now begin with puncturing the peripheral veins (fig. 6)

Emptying and cleaning the tubing system

- Close the flow control of the injection tube and disconnect the injection tube from the tubing system. Hold both ends into the cup.
- Take the infusion bottle with the remaining artificial blood from the stand and first remove the infusion system and then the rubber cork of the bottle.
- Open the flow control of the infusion system and let the artificial blood run off into the infusion bottle.
- Take the injection arm from the stand and hold it in an upright position so that the hand shows upward. Hold the end of the tube sealed with the white cap straight upwards so that the artificial blood can flow into the infusion bottle (fig. 7)
- Reseal the tube ends with the white caps and close the infusion bottle with the screw cap. The artificial blood can be reused and can also be further diluted with water, if necessary.

Exchanging the skin and the tubing system

- The core of the injection arm consists of 2 parts (arm core and hand core) that are connected by a guide device. When changing the skin and the tubing system, the hand core remains in the skin and only the arm core needs to be taken out.
- Remove the screw ring at the end of the injection arm, remove the white sealing caps from the tube ends and remove the white cover plates. Pull the tube end through the locating hole in the skin.
- Pull the skin off the arm core of the injection arm like a glove. Although the skin is very elastic, you should avoid extreme stretching.
- Unfasten the 5 Velcro strips holding the tubing system in place in the arm core and remove the tubing system.

English

Injection Arm

- Clean the arm core and the skin, rinsing off any residues of artificial blood with clear water.
- Place the new tubing system into the provided grooves and refasten it using the Velcro strips.
- Now hold the skin with the inner hand core pointing downwards and put some talcum powder into the opening. Keep the opening well sealed and spread the powder evenly by tilting in all directions. The powder improves the smoothness of the skin.
- Now you can pull the skin over the arm core again.
- Put the end of the tube back through the location hole in the skin and the white cover plates. The rim of the skin should be squeezed between the two cover plates. Now you can refasten everything with the screw ring and replace the white sealing caps onto the tube ends – finished!

Product numbers for replacement parts

XP102	Infusion bottle with suspension net
XP103	10 bags of artificial blood powder, 1 infusion system with air filter
XP104	3 complete tubing systems
XP106	Replacement skin with hand core
XP107	Arm core (without hand)

In case of complaints please indicate the serial number imprinted on the core of the injection arm under the white cover plates.

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